



The versatility you need

SureSigns VM6 portable bedside monitor

The Philips SureSigns VM6 now includes apnea alarming, adding to the numerous essential measurements already available in this easy-to-use, compact monitor.

When combined with additional measurements including user-selectable 3- or 5-lead ECG, pulse oximetry, non-invasive blood pressure (NBP), and continuous temperature, the SureSigns VM6 offers convenience, reliability, and quality patient care.

PHILIPS
sense and simplicity

SureSigns VM6. Beyond

Easy to learn, easy to use

The VM6, like all SureSigns monitors, is easy to set up, easy to learn, easy to use, and easy to support. An enhanced patient menu contains added patient ID input fields such as MRN, Transaction/Visit, Location and Operator. These IDs can easily be entered with a barcode scanner (optional).


Intuitive navigation puts your most commonly used features at your fingertips. And using the VM6 has become even easier with the addition of tool tips, import/export of configuration settings, and options for printing waveforms (7 or 20 seconds).

Connect to the SureSigns VSV: a simple way to be everywhere at once

Having a clear picture of each patient's status, at a single glance, can help relieve stress on care teams. The Philips SureSigns Vital Signs Viewer is an affordable, dependable, and simple way to add the benefits of central viewing and alarming without the complexity or expense of a traditional PC-based central station. Connect up to 12 SureSigns VM monitors at a time. It is one more way Philips is making it easier to free yourself to do exactly what you need to do to care for your patients.

All the support you need

The SureSigns VM6 is backed by high-quality technology and support. Software upgrades, via the USB port, are remarkably quick and easy. The sleek modular design simplifies set up and repair. Our warranty adds to your peace of mind.




Sleek roll stand with integrated basket and locking casters makes it easier to take the SureSigns VM6 patient monitor where it's needed. Wall mount also available. Each is designed for easy mounting and quick release.



The Philips VSV allows you to check the status of up to 12 patients at a glance.



the basics.



Integrated bedrail hanger for maximum flexibility.

Large 8.4" (20.3 cm) color display with up to four waveforms and numerics is bright and easy to read.

Sturdy handle and rugged housing for easy portability.

Lightweight at only 3.0 kg (6.6 lbs).

Lithium ion battery is standard.

Enhanced patient menu containing additional ID input fields.

Fixed keys for quick access to main functions.

Non-invasive blood pressure (NBP), including manual, interval, and STAT modes. Supports the Philips hardware module.

Integrated strip recorder for convenient documentation (optional).

Up to 96 hours of measurement trends.

Bar code scanner for convenient entry of patient ID (optional).

USB port for easy software upgrades.

User-selectable 3- and 5-lead ECG with pacemaker detection, including a "leads off" indicator.

Philips SpO₂ technology compatible with a wide variety of adult, pediatric, and neonatal sensors.*

* Adapter cable required.

Product specifications	
Width, height, depth	W: 26 cm (10.2"), H: 22 cm (8.6"), D: 14.5 cm (5.7")
Weight	3.0 kg (6.6 lbs), excluding optional recorder
Screen	8.4" (20.3 cm) SVGA TFT-AM LCD
NBP	<ul style="list-style-type: none"> Oscillometric using stepwise deflation pressure • Adult measurement range <ul style="list-style-type: none"> – Systolic: 30 – 270 mmHg – Diastolic: 10 – 245 mmHg – MAP: 20 – 255 mmHg – Pulse rate range: 40 – 300 bpm • Pediatric measurement range <ul style="list-style-type: none"> – Systolic: 30 – 180 mmHg – Diastolic: 10 – 150 mmHg – MAP: 20 – 160 mmHg – Pulse rate range: 40 – 300 bpm • Neonatal measurement range <ul style="list-style-type: none"> – Systolic: 30 – 130 mmHg – Diastolic: 10 – 100 mmHg – MAP: 20 – 120 mmHg – Pulse rate range: 40 – 300 bpm • Blood pressure accuracy <ul style="list-style-type: none"> – Maximum standard deviation: 8 mmHg – Maximum mean error: ± 5 mmHg • Pulse rate accuracy <ul style="list-style-type: none"> – 40 – 100 bpm: ± 5 bpm – 101 – 200 bpm: $\pm 5\%$ of reading – 201 – 300 bpm: $\pm 10\%$ of reading
IBP (optional)	<ul style="list-style-type: none"> • ABP/CVP/PAP measurement range: -40 – 360 mmHg • Overall accuracy (including transducer): greater of ± 4 mmHg or $\pm 4\%$
SpO ₂	<ul style="list-style-type: none"> • Measurement range: 0 – 100% • Accuracy depends on sensor
ECG	<ul style="list-style-type: none"> • 3- or 5-lead; heart rate range: 15 – 300 bpm • Heart rate accuracy: greater of $\pm 1\%$ or ± 5 bpm • Arrhythmia analysis
Impedance respiration	<ul style="list-style-type: none"> • Technique: Trans-thoracic impedance • Measurement range: 0 – 150 rpm • Accuracy: <ul style="list-style-type: none"> – ± 1 rpm in the range of 0 – 120 rpm – ± 2 rpm in the range of > 120 rpm
Temperature	<ul style="list-style-type: none"> • Measurement Range: 25 – 45°C (77 – 113°F) • Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
Battery	<ul style="list-style-type: none"> • 10.8 – 11.1V (lithium ion battery) • Operating time: 4 hours with continuous monitoring and NBPs every 15 minutes
Data output	• HL7 format, via Ethernet port
Patients	• Adult, pediatric, neonatal

863064	SureSigns VM6
863065	SureSigns VM6 with IBP
R01	Recorder (optional)
863067	SureSigns VSV

* Additional supplies available. Please contact your Philips representative for more information.
 ** Depending on model ordered.

Please visit www.philips.com/SureSigns



© 2011 Koninklijke Philips Electronics N.V.
 All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Supplies included*	To re-order
Adult antimicrobial NBP cuff	M4555A
NBP interconnect tubing	M1599B
Adult finger SpO ₂ sensor	M1196A
3-lead ECG lead set with grabber cable	989803143171 (IEC) or 989803143181 (AAMI)
Recorder paper (5 rolls)**	989803136891
Also available	To order
Roll stand with basket	989803144001
Wall mount	989803144011
Bar code scanner kit	1D: 989803167691 2D: 989803147821

A simple way to be everywhere at once

Connect with the SureSigns Vital Signs Viewer to show status from up to 12 patient monitors at a time.

Getting started

Each monitor comes with set-up guide, instructions for use, quick reference cards, user training DVD, and service manual on CD.



Supplies for the SureSigns VM6 include a wide selection of disposable and reusable SpO₂ sensors, all designed for accuracy, performance, and durability.



1800 US Hwy 51 N
Woodruff, WI 54568
715-356-0205

SpotMonitor.com
Monitors for Life
An Allied 100 Company

To order visit
SpotMonitor.com
or call 800-920-9919

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
 4522 962 70921 * APR 2011